







## GLOSHIELD® SELECTION AND INSTALLATION GUIDE

**Table 1:** GloShield® & Light Cable Compatibility

GloShield® Model	Light Cable Type	Light Cable Connector Profile
Storz® (Color: Sky Blue) Order: REF 16-3001 (Carton) Track: REF 16-2001 (Pouch) 	Storz® (495 Series and similar profile except 495 NCSC and NAC)  Stryker® (233-050-064, 233-050-084, and similar profile)  Conmed® (1050 w/ 7453 connector)	
Olympus® (Color: Orchid) Order: REF 17-3001 (Carton) Track: REF 17-2001 (Pouch) 	Olympus® (WA03200A, WA03210A, WA03310A, and similar profile)	



**Figure 1:** Proper positioning (left) and improper positioning (middle, right) of Sky Blue GloShield®



**Figure 2:** Proper positioning (left) and improper positioning (right) of Orchid GloShield®